

# Case study

## Protein A resin screening for Mouse IgG

### INTRODUCTION

A leading medical technology company, a pioneer in diagnostics, was looking for a process optimization to improve recovery and purification solution.

### CHALLENGE

The customer's purification process operated at reduced binding capacity to enable the purification of mouse IgG. To improve process efficiency, however, they hoped to identify a solution that would enable them to obtain equivalent purity with a higher dynamic binding capacity.

### AVANTOR SOLUTION

An Avantor team executed a study based on protocol developed collaboratively with the customer and using customer-provided resins, IgG and conditions. The study demonstrated the performance of PROchievA™ resin compared to that of other resins used for the purification of mouse IgG.

### RESULT – VALUE TO CUSTOMER

The results of the study indicated a significant increase in efficiency with the use of PROchievA™, with time saved in screening and development.

### CUSTOMER BENEFIT

Based on these results, the customer implemented a new process with PROchievA™ resin and achieved improved recovery and purity in their manufacturing process.



### Protein A resin screening for Mouse IgG

Protein A resin	Residence time (min)	Mouse IgG1 DBC (mg/mL)
PROchievA™	8	31.6
	6	30.3
	4	28.3
	2	20.3
Competitor resin 1	2	7.5
Competitor resin 2	2	3.7

### PROchievA™ ligand leaching data

The experiment was performed with a 1 ml column at 2 minute residence time using a pH 2.9 elution buffer, Leached protein A in eluate sample was quantified by ELISA.

Sample	Protein A concentration
Elution from 2 min residence time	1.21 ppm

### Outcome

Overall, PROchievA™ showed better performance, with higher purity and dynamic binding capacity.

